

# Getting Started with



## 1. Download

Download the CViConnect PRO app from your iPad's App Store.




## 2. Log In

Enter your unique username and password.




## 3. Select Student

Using the <For Student> drop down menu, select the student you wish to work with.




## 4. Select an Activity

Touch the play icon in the center of the Activity you'd like to open. Carefully read the instructions provided.




## 5. Preview Look Detection

Touch the more options icon (...), select **Preview Eye Detection**. Align the iPad's forward facing camera with your student's face.



## 6. Begin Activity

Press the play icon to begin the Activity.



## 7. Complete the Activity

To complete the Activity session, press the pause icon in the upper right-hand corner. Select **Complete** from the prompt. Lastly, in the text box provided, enter your environmental notes or other observations relevant to your student's progress.

